

**Amendments to the Claims:**

This listing of claims will replace all prior versions, and listings, of claims in the application:

1.       (*Previously Presented*) A medical device package comprising:  
a main cap member with a cavity therein, the main cap member including:  
a proximal end; and  
a distal end; and  
a minor cap member,  
wherein:  
the cavity has a cavity opening at the proximal end of the main cap member;  
the cavity is configured to receive, and to securely and removably retain, no more than a portion of a medical device with electrical contacts therein such that the electrical contacts project from the cavity opening and minor cap member, and  
the minor cap member is configured to seal the cavity opening once the medical device has been partially received in the cavity.
2.       (*Original*) The medical device package of claim 1, wherein the main cap member includes at least one lateral channel and wherein the medical device is securely and removably retained by a friction fit between the medical device and the at least one lateral channel.
3.       (*Original*) The medical device package of claim 1, wherein the main cap member includes a directional marker.
4.       Canceled
5.       Canceled
6.       (*Original*) The medical device package of claim 1, wherein the minor cap member is a breachable minor cap member.
7.       (*Original*) The medical device package of claim 1, wherein the minor cap member is configured for permanent attachment to the medical device.

8.        (*Original*) The medical device package of claim 1, wherein the main cap member further includes a distal end cavity configured for disabling the medical device.
9.        (*Previously Presented*) A medical device package kit comprising:  
a main cap member with a cavity therein, the main cap member including:  
a proximal end; and  
a distal end;  
a minor cap member; and  
a connector,  
wherein:  
the cavity has a cavity opening at the proximal end of the main cap member;  
the cavity is configured to receive, and to securely and removably retain, no more than a portion of a medical device with electrical contacts therein such that the electrical contacts project from the cavity opening and minor cap member;  
the minor cap member is configured to seal the cavity opening once the medical device has been partially received in the cavity; and  
the connector is configured to engage the medical device during removal of the medical device from the cavity.
10.       (*Original*) The medical device package kit of claim 9, wherein the connector includes a connector directional marker.
11.       (*Original*) The medical device package kit of claim 9, wherein the connector is further configured to breach the minor cap member.
12.       (*Previously Presented*) The medical device package kit of claim 9, wherein the connector includes at least one strip engaging element for contacting the electrical contacts of the medical device.
13.       (*Previously Presented*) A method for extracting a medical device from a medical device package, the method comprising:  
providing:  
a medical device package with a medical device partially contained therein, the medical device package including:

- a main cap member with a cavity therein, the main cap member including:
  - a proximal end; and
  - a distal end;
- a minor cap member; and

wherein:

- the cavity has a cavity opening at the proximal end of the main cap member;
- the cavity is configured to receive, and to securely and removably retain, no more than a portion of a medical device with electrical contacts therein such that the electrical contacts project from the cavity opening and minor cap member; and
- the minor cap member is sealing the cavity opening; and

a connector;

breaching the minor cap member with the connector such that at least a portion of the connector has entered the cavity of the medical device cavity;

engaging the medical device with the connector; and

extracting the connector and engaged medical device from the cavity of the medical device package.

14.     (*Original*)   The method of claim 13, wherein the medical device is an integrated medical device that includes a tissue penetration member and a test strip.

15.     (*Previously Presented*)   A method for extracting a medical device from a medical device package and subsequently disabling the medical device, the method comprising:

providing:

- a medical device package with a medical device contained therein at a first position, the medical device package including:

- a main cap member with a cavity therein, the main cap member including:
    - a proximal end; and
    - a distal end;
  - a minor cap member; and

wherein:

- the cavity has a cavity opening at the proximal end of the main cap member;
  - the cavity is configured to receive, and to securely and removably retain, no more than a portion of a medical device with electrical contacts therein such that the electrical contacts project from the cavity opening and minor cap member;

the minor cap member is sealing the cavity opening; and  
a connector;  
breaching the minor cap member with the connector such that at least a portion of the connector has entered the cavity of the medical device cavity;  
engaging the medical device with the connector;  
extracting the connector and engaged medical device from the cavity of the medical device package; and  
subsequently inserting the connector and engaged medical device into a cavity of the medical device package to a second position whereby the medical device is disabled from subsequent use.

16.     (*Original*) The method of claim 15, wherein the inserting step includes inserting the engaged medical device into the cavity of the main cap member.

17.     (*Original*) The method of claim 16, wherein the inserting step includes inserting the engaged medical device into a distal end cavity of the medical device package.

18.     (*Original*) The method of claim 16 further including:  
disengaging the connector from the disabled medical device; and  
withdrawing the connector from the medical device package.